

FICHA TÉCNICA

LIMAS K NI-TI PROCLINIC		
REFERENCIA	Ver anexo 1	
FABRICANTE / RESPONSABLE	SHENZHEN SUPERLINE TECHNOLOGY CO., LTD. Room 302, Nanshan Zhiyuan A4, No.1001 Xueyuan Avenue, Changyuan community, Taoyuan Street, Nanshan District, Shenzhen, 518000 Guangdong, P.R. China	
CLASIFICACIÓN DEL PRODUCTO	Producto Sanitario clase IIa, según la Directiva 93/42/CEE relativa a productos sanitarios.	
DESCRIPCIÓN	Dispositivos manuales de endodoncia para la práctica odontológica. No estériles.	
COMPOSICIÓN	Parte activa: Acero inoxidable 304H, Níquel-Titanio (NITI). Tope: Silicona. Mango: Tereftalato de polibutileno PBT.	
ENVASE	6 unidades.	
FINALIDAD PREVISTA	Dispositivos para canales radiculares endodónticos destinados a explorar, dar forma y limpiar los canales radiculares.	
LIMPIEZA Y ESTERILIZACIÓN	Los dispositivos se deben examinar, limpiar y esterilizar antes de cada uso. Usar detergentes y desinfectantes adecuados para instrumentos médicos o dentales. Esterilizar en autoclave. Desechar el producto si está dañado o disminuye su rendimiento.	
ALMACENAMIENTO Y TRANSPORTE	No requiere condiciones especiales de almacenamiento y transporte.	

TECHNICAL DATA SHEET

K NI-TI PROCLINIC FILES		
REFERENCES	See Annex 1.	
MANUFACTURER	SHENZHEN SUPERLINE TECHNOLOGY CO., LTD. Room 302, Nanshan Zhiyuan A4, No.1001 Xueyuan Avenue, Changyuan community, Taoyuan Street, Nanshan District, Shenzhen, 518000 Guangdong, P.R. China	
CLASSIFICATION	Class IIa medical device, according to Directive 93/42/EEC on medical devices.	
DESCRIPTION	Endodontic handpieces for dental practice. Non sterile.	
COMPOSITION	Active part: 304H stainless steel, Nickel-Titanium (NITI). Stop: Silicone. Handle: Polybutylene terephthalate PBT.	
PRESENTATION	6 units.	
INTENDED USE	Endodontic root canal devices for exploring, shaping and cleaning root canals.	
CLEANING AND STERILIZATION	Devices should be examined, cleaned and sterilized before each use. Use detergents and disinfectants suitable for medical or dental instruments. Sterilize in autoclave. Discard the product if it is damaged or its performance decreases.	
STORAGE AND TRANSPORT	Store in a cool and dry place. Keep away from sunlight. Protect from moisture and avoid excessive heat.	

ANEXO I

Ref.	Nombre / Name / Nome
17605	LIMAS NI-TI N.15 (21mm) PROCLINIC
17606	LIMAS NI-TI N.20 (21mm) PROCLINIC
17607	LIMAS NI-TI N.25 (21mm) PROCLINIC
17608	LIMAS NI-TI N.30 (21mm) PROCLINIC
17609	LIMAS NI-TI N.35 (21mm) PROCLINIC
17610	LIMAS NI-TI N.55 (21mm) PROCLINIC
17611	LIMAS NI-TI N.50 (21mm) PROCLINIC
17612	LIMAS NI-TI N.40 (21mm) PROCLINIC
17613	LIMAS NI-TI N.45 (21mm) PROCLINIC
17614	LIMAS NI-TI N.60 (21mm) PROCLINIC
17615	LIMAS NI-TI N.70 (21mm) PROCLINIC
17616	LIMAS NI-TI N.80 (21mm) PROCLINIC
17617	LIMAS NI-TI N.15 (25mm) PROCLINIC
17618	LIMAS NI-TI N.20 (25mm) PROCLINIC
17619	LIMAS NI-TI N.25 (25mm) PROCLINIC
17620	LIMAS NI-TI N.30 (25mm) PROCLINIC
17621	LIMAS NI-TI N.35 (25mm) PROCLINIC
17622	LIMAS NI-TI N.45 (25mm) PROCLINIC
17623	LIMAS NI-TI N.40 (25mm) PROCLINIC
17624	LIMAS NI-TI N.50 (25mm) PROCLINIC
17625	LIMAS NI-TI N.80 (25mm) PROCLINIC
17626	LIMAS NI-TI N.70 (25mm) PROCLINIC
17627	LIMAS NI-TI N.55 (25mm) PROCLINIC
17628	LIMAS NI-TI N.60 (25mm) PROCLINIC