

FICHA TÉCNICA

LÉNTULOS PROCLINIC		
REFERENCIAS	17775: LENTULOS N.25 PROCLINIC (25mm). 17776: LENTULOS N.30 PROCLINIC (25mm). 17777: LENTULOS N.35 PROCLINIC (25mm). 17778: LENTULOS N.40 PROCLINIC (25mm).	
FABRICANTE / RESPONSABLE	SHENZHEN SUPERLINE TECHNOLOGY CO., LTD. Room 302, Nanshan Zhiyuan A4, No.1001 Xueyuan Avenue, Changyuan community, Taoyuan Street, Nanshan District, Shenzhen, 518000 Guangdong, P.R. China	
CLASIFICACIÓN	Producto Sanitario Clase IIa, según la Directiva 93/42/CEE.	
DESCRIPCIÓN	-Instrumentos rotatorios para CA flexibles que facilitan el sellado y cementación de conductos radiculares. -Esterilizables en autoclave.	
ESPECIFICACIÓN	<u>Medidas:</u> 17775: LENTULOS N.25 PROCLINIC (25mm). 17776: LENTULOS N.30 PROCLINIC (25mm). 17777: LENTULOS N.35 PROCLINIC (25mm). 17778: LENTULOS N.40 PROCLINIC (25mm).	
COMPOSICIÓN	Acero inoxidable.	
ENVASE	Envase con 4 unidades.	
VIDA ÚTIL	La vida útil se estima en 3 años a partir de la fecha de fabricación, pero puede variar según las condiciones y frecuencia de uso y el número de ciclos de esterilización.	
USO PREVISTO	Facilitan la introducción del fondo de la cavidad, eliminando el aire existente. Sólo para uso por profesionales especializados en odontología.	
ALMACENAMIENTO Y TRANSPORTE	No requiere condiciones especiales de almacenamiento y transporte.	

TECHNICAL DATA SHEET

LENTULOS		
REFERENCES	17775: LENTULOS N.25 PROCLINIC (25mm). 17776: LENTULOS N.30 PROCLINIC (25mm). 17777: LENTULOS N.35 PROCLINIC (25mm). 17778: LENTULOS N.40 PROCLINIC (25mm).	
MANUFACTURER	SHENZHEN SUPERLINE TECHNOLOGY CO., LTD. Room 302, Nanshan Zhiyuan A4, No.1001 Xueyuan Avenue, Changyuan community, Taoyuan Street, Nanshan District, Shenzhen, 518000 Guangdong, P.R. China	
CLASSIFICATION	Class IIa Medical Device, according to Directive 93/42/EEC.	
DESCRIPTION	-Flexible rotary instruments for AC that facilitate the sealing and cementation of root canals. -Sterilizable in autoclave.	
SPECIFICATION	<u>Measures:</u> 17775: LENTULOS N.25 PROCLINIC (25mm). 17776: LENTULOS N.30 PROCLINIC (25mm). 17777: LENTULOS N.35 PROCLINIC (25mm). 17778: LENTULOS N.40 PROCLINIC (25mm).	
COMPOSITION	Stainless steel.	
PACKAGING	Package with 4 units.	
SHELF LIFE	The shelf life is estimated at 3 years from the date of manufacture but may vary depending on the conditions and frequency of use and the number of sterilization cycles.	
INTENDED USE	They facilitate the introduction of the bottom of the cavity, eliminating the existing air. These medical devices are used by professionals skilled in dentistry.	
STORAGE AND TRANSPORT	It does not require special storage and transportation conditions.	