

FICHA TÉCNICA

ADAPTADOR ASPIRASALIVA		
REFERENCIAS	82219: 6 a 11 mm. 82220: 11 a 16 mm.	
FABRICANTE / RESPONSABLE	TALLERES MESTRAITUA, S.L. Txorierrri Etorbidea, 60, 48150, Bizkaia ESPAÑA www.mestra.es	
DESCRIPCIÓN	Adaptador de eyectores para sistemas de aspiración de equipos dentales.	
ESPECIFICACIÓN	Adapta el aspirasaliva de 6 mm a la salida de la manguera de 11 mm. Adapta el aspirasaliva de 11 mm a la salida de la manguera de 16 mm. Características: Esterilizable en autoclave a 134°C.	
COMPOSICIÓN	Fabricado en material plástico.	
PRESENTACIÓN	82219: Envase con 5 unidades. 82220: Envase con 10 unidades.	
FINALIDAD PREVISTA	Conexión cánula aspirasaliva con manguera de aspiración.	
ALMACENAMIENTO Y TRANSPORTE	No requiere condiciones especiales de almacenamiento y transporte.	

TECHNICAL DATA SHEET

SALIVA EJECTOR ADAPTOR		
REFERENCES	82219: 6 a 11 mm. 82220: 11 a 16 mm.	
MANUFACTURER / RESPONSIBLE	TALLERES MESTRAITUA, S.L. Txorierri Etorbidea, 60, 48150, Bizkaia ESPAÑA www.mestra.es	
DESCRIPTION	Ejector adapter for dental equipment suction systems.	
SPECIFICATION	Adapts the 6 mm aspirasaliva to the 11 mm hose outlet. Adapts the 11 mm aspirasaliva to the 16 mm hose outlet. <u>Features:</u> Sterilizable in autoclave at 134°C.	
COMPOSITION	Made of plastic material.	
PACKAGE	82219: Package with 5 units. 82220: Package with 10 units.	
INTENDED USE	Suction cannula connection with suction hose.	
STORAGE AND TRANSPORT	It does not require special storage and transportation conditions.	