


PROCESS VALIDATION PRINTERS

CERTIFIED
VALIDATED
RELIABLE



LATEST
VERSION
WWW.DETAX.COM



	385 nm													405 nm										
Qualification	ASIGA Max / Mini	ASIGA Picc2	ASIGA PRO2	ASIGA PRO 4K	MICROLAY Versus	Miicraft Prime / Hyper Series	Miicraft Ultra Series	Miicraft Profession / Advance Series	Rapid shape D10 / D20 Series	Rapid shape D30 / D40 Series	Rapid shape D10 / D30 Series	Rapid shape D100	Strumann P series	W2P	Ackuretta	Flashforge Hunter	Microsay Eye Pro	Moonray S100	Phrozen Sonic 4K / XL	Prusa Research MEDICAL ONE	Shining 3D Accu-Fab D1s	Shining 3D Accu-Fab L4D	Sprintray Pro S	
temp	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	●	✓	✓	●	□	□	□	□	□	□	□	□	□
crown	✓	●	●	✓	✓	□	□	□	✓	✓	□	●	✓	●	●	●	●	●	●	●	●	●	●	●
denture	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	□	✓	✓	●	□	□	□	□	□	□	□	□	□
ortho	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	●	✓	✓	●	□	●	●	●	●	●	●	●	●
splint 2.0	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	●	✓	✓	●	●	●	●	●	●	●	●	●	●
splintmaster	□	□	□	□	□	□	□	□	□	□	□	●	□	●	●	●	●	●	●	●	●	●	●	●
IBT	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	●	✓	✓	●	●	●	●	●	●	●	●	●	●
tray 2.0	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	●	✓	✓	●	●	□	●	□	□	□	□	□	●
tryin	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	●	✓	✓	●	●	●	●	●	□	□	□	□	●

model	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	✓	✓	✓	✓	✓	✓	✓	✓	✓
model 2.0	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	□	✓	✓	□	✓	✓	✓	✓	□	□	□	□	✓
model T	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	□	✓	✓	□	✓	✓	✓	✓	□	✓	✓	✓	✓
model WW	✓	✓	✓	✓	□	□	□	□	□	□	□	□	●	□	✓	✓	✓	✓	✓	●	●	●	●	✓
gingiva	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	●	✓	✓	□	✓	✓	✓	✓	□	✓	✓	✓	✓
cast 2.0	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	●	✓	✓	□	✓	✓	✓	✓	□	✓	✓	✓	✓
Primeprint (made by DETAX)																								✓

edition: 03.02.23

FREEPRINT®
MED

FREEPRINT®
TEC

PRIMEPRINT
MED & TEC

PROCESS VALIDATION CURING DEVICES

GOOD TO KNOW ...


Qualification	LED														Bentley Strona Primeprint PPU	
	KK Optik Otoflash Q171 N2	Otoflash 250	Otoflash 500	Rapid Shape RS cure	DentalFarm Photop	Mecatronicore ES-Cure	Scheu Impress Cure	Phrozen Phrozen Cure	Sprint Ray Pro Cure	Formlabs Form Cure	Hey Gears PCU 3.0	Drove PCU LED N2	Shining 3D FacCure	Prusa Research Medical CW one		
temp	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
crown	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
denture	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
ortho	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
splint 2.0	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
splintmaster Q4/23	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
IBT	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
tray 2.0	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
tryin Q3/23	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
model	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
model 2.0	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
model T	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
model WW	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
gingiva	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
cast 2.0	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
Primeprint (made by DETAx)	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>

edition: 03.02.23



BOTTLE ROLLER

By using a bottle roller, optimum mixing of the material is achieved, thus preventing possible segregation. The Eco Bags can be homogenized with an appropriate attachment.




CLEANING

Best cleaning results of the production jobs are achieved when the pre- and post-cleaning are carried out in separate tanks in the ultrasonic unit. After cleaning with isopropanol, it is recommended to clean the bores/openings with compressed air.



POST-CURING UNIT

The post-curing units recommended in the instructions for use ensure optimum through-hardening and surface curing, thus a biocompatible end product, and ensure high color brilliance and transparency, without discoloration.



DETAx EXPERTS@

