

GC Initial™ REVO MC



This chart is intended as a reference guide only.



GC Initial™ REVO MC



16 A-D Keyshades	A					B					C				D		
	Opaque	Light	Light	Medium	Medium	Dark	Light	Light	Medium	Dark	Light	Light	Medium	Dark	Light	Medium	Dark
Body Light	100%	75%	50%	25%	0%	100%	75%	50%	0%	100%	75%	50%	0%	100%	50%	0%	
Body Dark	0%	25%	50%	75%	100%	0%	25%	50%	100%	0%	25%	50%	100%	0%	50%	100%	
Enamel	Light	Light	Medium	Medium	Dark	Light	Medium	Medium	Medium	Light	Medium	Medium	Dark	Medium	Medium	Medium	

BLEACH	
Opaque	Bleach
Body Bleach	100%
Enamel	Bleach

Optimal & individual usage	SHOULDER																	
	Body Mixture	50%	50%	50%	50%	50%	50%	50%	50%	50%	50%	50%	50%	50%	50%	50%	50%	50%
	Shoulder Universal	50%	50%	50%	50%	50%	50%	50%	50%	50%	50%	50%	50%	50%	50%	50%	50%	
	EFFECT																	
Clear Fluorescent	For optimal and individual usage.																	
Opal Booster																		
Translucent Neutral																		

SHOULDER	
Body Bleach	50%
Shoulder Universal	50%

EFFECT	
Clear Fluorescent	For optimal and individual usage.
Opal Booster	
Translucent Neutral	

GUM	
Opaque	Gum Universal
Gum Light	Individual usage
Gum Dark	Individual usage

