

KAVO DENTAL HOSES

The dental hoses for Tecnomed Italia dental units are made of silicone; they are Latex Free and comply with the UNI CEI EN ISO 13485 quality system.

The chemically stable gold-plated contacts ensure high electrical conductivity and maximum handpiece performance.

The supply tubes comply with REGULATION (EU) No. 10/2011 concerning plastic materials objects intended to meet food products and the exhaust tubes are resistant to the lubrication oils.

Tecnomed Italia dental hoses meet the strictest medical standards, complying with the essential requirements of the European Directives (93/42 / EEC) 2007/47 / EC on medical devices and 2000/95 / EC RoHs on the restriction of the use of certain specific substances in electrical and electronic equipment.

A complete range of dental hoses:

- Turbines
- Micromotors
- Scalers
- Syringes
- Freely Adaptable Dental Hoses

Features

- ⊕ High quality silicone sheath.
- ⊕ Latex free.
- ⊕ 12 months warranty
















Technical data

- ⊕ Maximum pressure supply/return air 4,5 bar / 65 psi
- ⊕ Maximum pressure air/water spray 4 bar / 58 psi
- ⊕ Tensile strenght 70N

MICROMOTORS TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR K700 K701 K703	1057 1060		<p>JUMPER 3-4</p> <p>Spray regulation Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10,5 mm</p> <p>1350 mm</p>		PC148KN13
MICROMOTOR K700 K701 K703	1058 1062 1063 1065 1066 1078 1080		<p>JUMPER 1-3-4</p> <p>Spray regulation Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10,5 mm</p> <p>1350 mm</p>		PC148KN
MICROMOTOR K700 K701 K703	1058 1062 1063 1065 1066 1078 1080		<p>JUMPER 1-3-4</p> <p>Spray regulation Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1350 mm - Ø 10,5 mm</p> <p>1450 mm</p>		PC148KN14










MICROMOTORS TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR K700 K701 K703	E30S - E30T E50S - E70S E50TM - E70C E70T - E80C E80T 1058 LIFE		<p>RES. 4,2 KΩ</p> <p>Spray regulation  Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10,5 mm</p> <p>1350 mm</p> 		PC148KN5
MICROMOTOR K700 K701 K703	E30S E30T		<p>RES. 4,2 KΩ</p> <p>Spray regulation  Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10,5 mm</p> <p>1350 mm</p> 		PC148KN3
MICROMOTOR K190 K196 K200 TMI MJ660	1057 1060		<p>JUMPER 3-4</p> <p>Spray regulation  Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10,5 mm</p> <p>1350 mm</p> 		PC146KN




MICROMOTORS TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR K190 K196 K200 TMI MJ660	1058 1062 1063 1065 1066 1080		<p>JUMPER 1-3-4 Spray regulation Swivel LIGHT GRAY RAL7035 1250 mm - Ø 10,5 mm 1350 mm</p>		PC135KN
MICROMOTOR K190 K196 K200 TMI MJ660	1066 1080		<p>JUMPER 1-3-4 Spray regulation Swivel LIGHT GRAY RAL7035 1400 mm - Ø 10,5 mm 1500 mm</p>		PC135KNL










HIGHSPEED TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
HIGHSPEED HSH	1057 1058 1060 1062 1063 1065 1066 1078 1080		 <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10 mm</p> <p>1350 mm</p>		PC103KN
HIGHSPEED HSH	E50TM E70C E70T E80C E80T 1058 LIFE		 <p>RES. 2,2 KΩ RID SENSOR</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10 mm</p> <p>1350 mm</p>		PC104KN
HIGHSPEED HSH	E50T E70T E80T PENDED		 <p>RES. 2,2 KΩ RID SENSOR</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1550 mm - Ø 10 mm</p> <p>1650 mm</p>		PC104KN14










HIGHSPEED TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
<p>HIGHSPEED HSH</p>	<p>E30S E30T E50S E70S 1058 LIFE S</p>		<p>RES. 2,2 KΩ</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10 mm</p> <p>1350 mm</p> 		<p>PC105KN</p>

SYRINGE TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SYRINGE 3F KAVO (07730500)	1058 1062 1063 1065 1066 1080		 <p>LIGHT GRAY RAL7035</p> <p>1150 mm - Ø 10,5 mm</p> <p>1250 mm</p>		PC168KE
SYRINGE 4F KAVO (07730000)	1058 1062 1063 1065 1066 1080		 <p>LIGHT GRAY RAL7035</p> <p>1150 mm - Ø 10,5 mm</p> <p>1250 mm</p>		PC168KN
SYRINGE 4F KAVO (07730000)	1066 1080		 <p>LIGHT GRAY RAL7035</p> <p>1350 mm - Ø 10,5 mm</p> <p>1450 mm</p>		PC168KNL







SYRINGE TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGTH / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
<p>SYRINGE</p> <p>4F KAVO (07730000)</p>	<p>E50TM E70C E70T E80C E80T 1058 LIFE TM</p>		<p>RES. 750 Ω</p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10,5 mm</p> <p>1350 mm</p> 		PC168KV
<p>SYRINGE</p> <p>4F KAVO (07730000)</p>	<p>E50S E70S E70T E80T</p>		<p>RES. 750 Ω</p> <p>LIGHT GRAY RAL7035</p> <p>1350 mm - Ø 10,5 mm</p> <p>1450 mm</p> 		PC168KVL
<p>SYRINGE</p> <p>6F LUZZANI O.F. MINIBRIGHT INOX</p>	<p>E50S E70S E70T E80T</p>		<p>RES. 750 Ω</p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10,5 mm</p> <p>1495 mm</p> 		PC168KMB





SYRINGE TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SYRINGE LUZZANI 3F MINILIGHT INOX	1057 1058 1060 1062 1063 1065 1066 1078 1080		 1200 mm - Ø 10,5 mm 1445 mm LIGHT GRAY RAL7035		PC179KN
SYRINGE LUZZANI 3F MINILIGHT INOX	E50TM E70C E70T E80C E80T 1058 LIFE TM		 1250 mm - Ø 10,5 mm 1495 mm LIGHT GRAY RAL7035 RES. 750 Ω		PC179KN2
SYRINGE LUZZANI 3F MINILIGHT INOX	E50S E70S E70T E80T		 1500 mm - Ø 10,5 mm 1745 mm LIGHT GRAY RAL7035 RES. 750 Ω		PC179KN2L

SCALER TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SCALER* PIEZOLED EMS O. F.	E50S E50TM E70C E70S E70T E80C E80T 1058 LIFE S 1058 LIFE TM 1058LIFE CART		<p>RES. 1,3 KΩ RID SENSOR</p> <p>Spray regulation  Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 8 mm</p> <p>1350 mm</p>		PC151KRWG
SCALER EMS O. F.	E30S E30T		<p>RES. 1,3 KΩ</p> <p>Spray regulation  Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 8 mm</p> <p>1350 mm</p>		PC151KE
SCALER EMS O. F.	E30S E30T		<p>RES. 1,3 KΩ</p> <p>Spray regulation  Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1440 mm - Ø 8 mm</p> <p>1520 mm</p>		PC151KEL

SCALER TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SCALER EMS O. F.	E50S E50TM E70C E70S E70T E80C E80T 1058 LIFE S 1058 LIFE TM 1058LIFE CART		<p>RES. 1,3 KΩ RID SENSOR</p> <p>Spray regulation  </p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 8 mm</p> <p>1350 mm</p>		PC151KEW
SCALER EMS O. F.	1058 1062 1063 1065 1066 1080		<p>JUMPER 1-3</p> <p>Spray regulation  </p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 8 mm</p> <p>1350 mm</p>		PC151KEWG

SCALER TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SCALER * EMS O. F.	E50S E50TM E70C E70S			<p>REVERSE POLARITY</p>	PC152KEWG
					PC152KEWGI
KAVO PIEZOLED	E70T E80C E80T 1058				
SCALER * EMS O. F.	E30S E30T				PC153KRWG
KAVO PIEZOLED					

