

KAVO DENTAL HOSES

The dental hoses for Tecnomed Italia dental units are made of silicone; they are Latex Free and comply with the UNI CEI EN ISO 13485 quality system.

The chemically stable gold-plated contacts ensure high electrical conductivity and maximum handpiece performance.

The supply tubes comply with REGULATION (EU) No. 10/2011 concerning plastic materials objects intended to meet food products and the exhaust tubes are resistant to the lubrication oils.

Tecnomed Italia dental hoses meet the strictest medical standards, complying with the essential requirements of the European Directives (93/42 / EEC) 2007/47 / EC on medical devices and 2000/95 / EC RoHs on the restriction of the use of certain specific substances in electrical and electronic equipment.

A complete range of dental hoses:

- Turbines
- Micromotors
- Scalers
- Syringes
- Freely Adaptable Dental Hoses

Features

- ⊕ High quality silicone sheath.
- ⊕ Latex free.
- ⊕ 12 months warranty

Technical data

- ⊕ Maximum pressure supply/return air 4,5 bar / 65 psi
- ⊕ Maximum pressure air/water spray 4 bar / 58 psi
- ⊕ Tensile strenght 70N



MICROMOTORS TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR K700 K701 K703	1057 1060		JUMPER 3-4 Spray regulation Swivel LIGHT GRAY RAL7035 1250 mm - Ø 10,5 mm 1350 mm		PC148KN13
MICROMOTOR K700 K701 K703	1058 1062 1063 1065 1066 1078 1080		JUMPER 1-3-4 Spray regulation Swivel LIGHT GRAY RAL7035 1250 mm - Ø 10,5 mm 1350 mm		PC148KN
MICROMOTOR K700 K701 K703	1058 1062 1063 1065 1066 1078 1080		JUMPER 1-3-4 Spray regulation Swivel LIGHT GRAY RAL7035 1350 mm - Ø 10,5 mm 1450 mm		PC148KN14

MICROMOTORS TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR K700 K701 K703	E30S - E30T E50S - E70S E50TM - E70C E70T - E80C E80T 1058 LIFE		RES. 4,2 KΩ Spray regulation  Swivel  LIGHT GRAY RAL7035 1250 mm - Ø 10,5 mm 1350 mm		PC148KN5
MICROMOTOR K700 K701 K703	E30S E30T		RES. 4,2 KΩ Spray regulation  Swivel  LIGHT GRAY RAL7035 1250 mm - Ø 10,5 mm 1350 mm		PC148KN3
MICROMOTOR K190 K196 K200 TMI MJ660	1057 1060		JUMPER 3-4 Spray regulation  Swivel  LIGHT GRAY RAL7035 1250 mm - Ø 10,5 mm 1350 mm		PC146KN

MICROMOTORS TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR K190 K196 K200 TMI MJ660	1058 1062 1063 1065 1066 1080				PC135KN
MICROMOTOR K190 K196 K200 TMI MJ660	1066 1080				PC135KNL







HIGHSPEED TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGTH / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
HIGHSPEED HSH	1057		Spray regulation Swivel LIGHT GRAY RAL7035 1250 mm - Ø 10 mm 1350 mm		PC103KN
	1058				
	1060				
	1062				
	1063				
	1065				
	1066				
	1078				
	1080				
HIGHSPEED HSH	E50TM		RES. 2,2 KΩ RID SENSOR Spray regulation Swivel LIGHT GRAY RAL7035 1250 mm - Ø 10 mm 1350 mm		PC104KN
	E70C				
	E70T				
	E80C				
	E80T				
	1058 LIFE				
HIGHSPEED HSH	E50T		RES. 2,2 KΩ RID SENSOR Spray regulation Swivel LIGHT GRAY RAL7035 1550 mm - Ø 10 mm 1650 mm		PC104KN14
	E70T				
	E80T PENDED				






HIGHSPEED TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
HIGHSPEED HSH	E30S E30T E50S E70S 1058 LIFE S		RES. 2,2 KΩ Spray regulation Swivel LIGHT GRAY RAL7035 1250 mm - Ø 10 mm 1350 mm		PC105KN




SYRINGE TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SYRINGE 3F KAVO (07730500)	1058 1062 1063 1065 1066 1080				PC168KE
SYRINGE 4F KAVO (07730000)	1058 1062 1063 1065 1066 1080				PC168KN
SYRINGE 4F KAVO (07730000)	1066 1080				PC168KNL

SYRINGE TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGTH / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SYRINGE 4F KAVO (07730000)	E50TM E70C E70T E80C E80T 1058 LIFE TM		RES. 750 Ω 		PC168KV
SYRINGE 4F KAVO (07730000)	E50S E70S E70T E80T		RES. 750 Ω 		PC168KVL
SYRINGE 6F LUZZANI O.F. MINIBRIGHT INOX	E50S E70S E70T E80T		RES. 750 Ω 		PC168KMB



SYRINGE TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SYRINGE LUZZANI 3F MINILIGHT INOX	1057 1058 1060 1062 1063 1065 1066 1078 1080				PC179KN
SYRINGE LUZZANI 3F MINILIGHT INOX	E50TM E70C E70T E80C E80T 1058 LIFE TM				PC179KN2
SYRINGE LUZZANI 3F MINILIGHT INOX	E50S E70S E70T E80T				PC179KN2L

SCALER TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGTH / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SCALER PIEZOLED	E50S E50TM E70C E70S E70T E80C E80T 1058 LIFE S 1058 LIFE TM 1058LIFE CART		RES. 1,3 KΩ RID SENSOR Spray regulation Swivel LIGHT GRAY RAL7035 1250 mm - $\varnothing$ 8 mm 1350 mm		PC151KRWG
SCALER EMS O. F.	E30S E30T		RES. 1,3 KΩ Spray regulation Swivel LIGHT GRAY RAL7035 1250 mm - $\varnothing$ 8 mm 1350 mm		PC151KE
SCALER EMS O. F.	E30S E30T		RES. 1,3 KΩ Spray regulation Swivel LIGHT GRAY RAL7035 1440 mm - $\varnothing$ 8 mm 1520 mm		PC151KEL

SCALER TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGTH / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SCALER EMS O. F.	E30S E30T E50S		RES. 1,3 KΩ Spray regulation  Swivel  	 	PC151KEW
SCALER EMS O. F.	1058 1062 1063 1065 1066 1080		JUMPER 1-3 Spray regulation  Swivel  	 	PC151KEWG