



FICHA TÉCNICA

PORTA ALICATES ALUMINIO Y PLASTICO MARCA: PROCLINIC		
REFERENCIAS	L6032	
FABRICANTE / RESPONSABLE	G&H Wire Company 2165 Earlywood Drive, Franklin IN, 46131 USA Emergency Contact: 317-346-6655	
DESCRIPCIÓN	Porta alicates, para aproximadamente 20 instrumentos.	
COMPOSICIÓN	Barras de aluminio con extremos de plástico.	
PRESENTACIÓN	1 porta alicates por bolsa.	
FINALIDAD	Diseñado especialmente para alicates de ortodoncia.	
INSTRUCCIONES DE USO	Ver instrucciones en el producto.	
ALMACENAMIENTO Y TRANSPORTE	No se requieren condiciones especiales de almacenamiento y transporte.	

FICHA TÉCNICA

ALUMINUM AND PLASTIC PLIERS HOLDER BRAND: PROCLINIC		
REFERENCES	L6032	
MANUFACTURER	G&H Wire Company 2165 Earlywood Drive, Franklin IN, 46131 USA Emergency Contact: 317-346-6655	
DESCRIPTION	Pliers holder, for approximately 20 instruments.	
COMPOSITION	Aluminum bars with plastic ends.	
PACKAGING	1 organizer per bag.	
INTENDED USE	Specially designed for orthodontic pliers.	
INSTRUCTIONS FOR USE	See instructions on the product.	
HANDLING AND STORAGE CONDITIONS	No special storage and transport conditions are required.	