




RS RADIX-STIFTE



C 0100  **RS RADIX-STIFTE**



**Symbol**

  →  
p.73

Grade 5-DIN 3.7165%,  
Al 5.5-6.75%,  
V 3.5-4.5%,  
Fe 0.4%, C 0.2%,  
H 0.015%, Co 0.1%,  
N 0.05%, Ti rest




C 0101 **PENETRATION DRILL**  
ELARGISSEUR – KANALERWEITERER – FRESA DE PENETRACIÓN – FRESA DI PENETRAZIONE

N°	ISO	SIZE	Symbol	
x6	1	2	3	




C 0102 **PRECISION SPIRAL DRILL**  
FORET DE PRÉCISION – PŘÁZISONS-SPIRALBOHRER – FRESA DE PRECISIÓN – FRESA DI PRECISIONE

N°	ISO	SIZE	Symbol	
x6	1	2	3	




C 0103 **GAUGE**  
JAUGE – BOHRLEHRE – CALIBRATOR – CALIBRO


N°	ISO	SIZE	Symbol	
x1	1	2	3	



C 0104 **WRENCH**  
CLÉ – SCHLÜSSEL – LLAVE – CACCIAVITE

N°	ISO	SIZE	Symbol	
x1	1	2	3	

C 0099    →



C 0098    →



RS RADIX-STIFTE TRAY

C 105P 1x 