

# OPEN END MICROMOTOR DENTAL HOSES

The dental hoses for Tecnomed Italia dental units are made of silicone; they are Latex Free and comply with the UNI CEI EN ISO 13485 quality system.

The chemically stable gold-plated contacts ensure high electrical conductivity and maximum handpiece performance.

The supply tubes comply with REGULATION (EU) No. 10/2011 concerning plastic materials objects intended to meet food products and the exhaust tubes are resistant to the lubrication oils.

Tecnomed Italia dental hoses meet the strictest medical standards, complying with the essential requirements of the European Directives (93/42 / EEC) 2007/47 / EC on medical devices and 2000/95 / EC RoHs on the restriction of the use of certain specific substances in electrical and electronic equipment.

A complete range of dental hoses:

- Turbines
- Micromotors
- Scalers
- Syringes
- Freely Adaptable Dental Hoses










## Features

- ⊕ High quality silicone sheath.
- ⊕ Latex free.
- ⊕ 12 months warranty










## Technical data

- ⊕ Maximum pressure supply/return air 4,5 bar / 65 psi
- ⊕ Maximum pressure air/water spray 4 bar / 58 psi
- ⊕ Tensile strenght 70N




# BIEN AIR MICROMOTORS

| HANDPIECE TYPE                | DENTAL UNIT MODEL | IMAGE  | LENGHT / DIAMETER / COLOR   | CONNECTIONS   | ORDER CODE |
|-------------------------------|-------------------|--|---|---|------------|
| MICROMOTOR<br>MC2<br>(NO LUX) |                   |    |  <p>Spray regulation Swivel</p> <p>EXTENSION STRING LIGHT GRAY RAL7035</p> <p>1500 mm - Ø 10 mm</p>    |    | PC722      |
| MICROMOTOR<br>MC2             |                   |   |  <p>Spray regulation Swivel</p> <p>EXTENSION STRING LIGHT GRAY RAL7035</p> <p>1500 mm - Ø 10 mm</p>   |   | PC732      |
| MICROMOTOR<br>MC3             |                   |  |  <p>Spray regulation Swivel</p> <p>EXTENSION STRING ITAL GRAY RAL7038</p> <p>1650 mm - Ø 10,5 mm</p> |  | PC666      |

# BIEN AIR MICROMOTORS

| HANDPIECE TYPE           | DENTAL UNIT MODEL | IMAGE  | LENGHT / DIAMETER / COLOR  | CONNECTIONS   | ORDER CODE |
|--------------------------|-------------------|--|--|---|------------|
| MICROMOTOR<br>MC3<br>HSH |                   |    |  <p>Spray regulation Swivel</p> <p>EXTENSION STRING</p> <p>LIGHT GRAY RAL7035</p> <p>1500 mm - Ø 11 mm</p>                    |    | PC666TM    |
| MICROMOTOR<br>MC3        |                   |   |  <p>Spray regulation Swivel Safety Handpiece</p> <p>EXTENSION STRING</p> <p>ITAL GRAY RAL7038</p> <p>1650 mm - Ø 10,5 mm</p> |   | PC666SH    |
| MICROMOTOR<br>MX2        |                   |  |  <p>Spray regulation Swivel</p> <p>EXTENSION STRING</p> <p>ITAL GRAY RAL7038</p> <p>1500 mm - Ø 9,5 mm</p>                  |  | PC734BXG   |

# BIEN AIR MICROMOTORS

| HANDPIECE<br>TYPE     | DENTAL UNIT<br>MODEL | IMAGE  | LENGHT / DIAMETER / COLOR  | CONNECTIONS   | ORDER<br>CODE |
|-----------------------|----------------------|--|--|---|---------------|
| MICROMOTOR<br><br>MCX |                      |  |  <p>Spray regulation</p> <p>Swivel</p> <p>EXTENSION STRING</p> <p>LIGHT GRAY RAL7038</p> <p>1750 mm - Ø 10,5 mm</p> |  | PC735MCX      |










# KAVO MICROMOTORS

| HANDPIECE TYPE                     | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR  | CONNECTIONS | ORDER CODE |
|------------------------------------|-------------------|-------|--|-------------|------------|
| MICROMOTOR<br>K192                 |                   |       | <p>EXTENSION STRING</p> <p>LIGHT GRAY RAL7035</p> <p>2000 mm - Ø 10,5 mm</p> |             | PC124PKN   |
| MICROMOTOR<br>K190<br>K196<br>K201 |                   |       | <p>EXTENSION STRING</p> <p>LIGHT GRAY RAL7035</p> <p>1600 mm - Ø 10,5 mm</p> |             | PC135PKN   |
| MICROMOTOR<br>K700<br>K701<br>K703 |                   |       | <p>EXTENSION STRING</p> <p>LIGHT GRAY RAL7035</p> <p>1600 mm - Ø 10,5 mm</p> |             | PC148PKN   |




# NSK MICROMOTORS

| HANDPIECE TYPE                                    | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR  | CONNECTIONS | ORDER CODE |
|---|-------------------|-------|--|-------------|------------|
| MICROMOTOR<br>NLX NANO - PLUS<br>ADEC<br>AE53-NLZ |                   |       | <p>LIGHT GRAY RAL7035</p> <p>1800 mm - Ø 8,5 mm</p> <p>Spray regulation Swivel</p>                         |             | PC601N     |
| MICROMOTOR<br>TIM 40<br>EL400                     |                   |       | <p>EXTENSION STRING</p> <p>LIGHT GRAY RAL7035</p> <p>1500 mm - Ø 9,5 mm</p> <p>Spray regulation Swivel</p> |             | PC602N     |
| MICROMOTOR<br>M40<br>HSH                          |                   |       | <p>EXTENSION STRING</p> <p>LIGHT GRAY RAL7035</p> <p>1600 mm - Ø 11 mm</p> <p>Spray regulation Swivel</p>  |             | PC666TM    |



# NSK MICROMOTORS

| HANDPIECE TYPE           | DENTAL UNIT MODEL | IMAGE  | LENGHT / DIAMETER / COLOR  | CONNECTIONS   | ORDER CODE |
|--------------------------|-------------------|--|--|---|------------|
| MICROMOTOR<br>M40<br>HSH |                   |    |  <p>Spray regulation Swivel Safety Handpiece</p> <p>EXTENSION STRING ITAL GRAY RAL7038</p> <p>1700 mm - Ø 11 mm</p> |    | PC901NSH   |
| MICROMOTOR<br>NLX-NL400  |                   |   |  <p>Spray regulation Swivel</p> <p>EXTENSION STRING LIGHT GRAY RAL7035</p> <p>1600 mm - Ø 10,5 mm</p>              |   | PC603N     |
| MICROMOTOR<br>NBX        |                   |  |  <p>Spray regulation Swivel</p> <p>EXTENSION STRING LIGHT GRAY RAL7035</p> <p>2200 mm - Ø 9 mm</p>                |  | PC604N     |







# FONA MICROMOTORS

| HANDPIECE<br>TYPE    | DENTAL UNIT<br>MODEL | IMAGE  | LENGHT / DIAMETER / COLOR  | CONNECTIONS   | ORDER<br>CODE |
|----------------------|----------------------|--|--|---|---------------|
| MICROMOTOR<br><br>BL |                      |  |  <p>Spray regulation</p> <p>Swivel</p> <p>EXTENSION STRING</p> <p>LIGHT GRAY RAL 7035</p> <p>1600 mm - Ø 9 mm</p> |  | PC233F        |

# PLANMECA MICROMOTORS

| HANDPIECE<br>TYPE     | DENTAL UNIT<br>MODEL | IMAGE  | LENGHT / DIAMETER / COLOR   | CONNECTIONS   | ORDER<br>CODE |
|-----------------------|----------------------|--|---|---|---------------|
| MICROMOTOR<br>MINETTO |                      |  | <p>RES. 2,87 KΩ</p> <p>Spray regulation</p> <p>Swivel</p> <p><a href="https://vimeo.com/480812578">https://vimeo.com/480812578</a></p> <p>EXTENSION STRING</p> <p>LIGHT GRAY RAL7035</p> <p>1800 mm - Ø 10 mm</p> |  | PC888P        |

# TEKNEDENTAL MICROMOTORS

| HANDPIECE TYPE                         | DENTAL UNIT MODEL | IMAGE   | LENGHT / DIAMETER / COLOR   | CONNECTIONS  | ORDER CODE |
|--|-------------------|---|---|--|------------|
| MICROMOTOR<br>TKD<br>DEFINITIVE<br>LED |                   |   |  <p>Spray regulation Swivel</p> <p>EXTENSION STRING ITAL GRAY RAL7038</p> <p>1650 mm - Ø 11 mm</p>                   |   | PC668      |
| MICROMOTOR<br>TKD<br>DEFINITIVE<br>LED |                   |  |  <p>Spray regulation Swivel Safety Handpiece</p> <p>EXTENSION STRING ITAL GRAY RAL7038</p> <p>1650 mm - Ø 11 mm</p> |  | PC668SH    |