

FREE MICROMOTOR DENTAL HOSES

The dental hoses for Tecnomed Italia dental units are made of silicone; they are Latex Free and comply with the UNI CEI EN ISO 13485 quality system.

The chemically stable gold-plated contacts ensure high electrical conductivity and maximum handpiece performance.

The supply tubes comply with REGULATION (EU) No. 10/2011 concerning plastic materials objects intended to meet food products and the exhaust tubes are resistant to the lubrication oils.

Tecnomed Italia dental hoses meet the strictest medical standards, complying with the essential requirements of the European Directives (93/42 / EEC) 2007/47 / EC on medical devices and 2000/95 / EC RoHs on the restriction of the use of certain specific substances in electrical and electronic equipment.

A complete range of dental hoses:

- Turbines
- Micromotors
- Scalers
- Syringes
- Freely Adaptable Dental Hoses








Features

- ⊕ High quality silicone sheath.
- ⊕ Latex free.
- ⊕ 12 months warranty







Technical data

- ⊕ Maximum pressure supply/return air 4,5 bar / 65 psi
- ⊕ Maximum pressure air/water spray 4 bar / 58 psi
- ⊕ Tensile strenght 70N







BIENAIR MICROMOTORS

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR MC2 (NO LUX)			 Spray regulation Swivel EXTENSION STRING LIGHT GRAY RAL7035 1500 mm - Ø 10 mm		PC722
MICROMOTOR MC2			 Spray regulation Swivel EXTENSION STRING LIGHT GRAY RAL7035 1500 mm - Ø 10 mm		PC732
MICROMOTOR MC3			 Spray regulation Swivel EXTENSION STRING ITAL GRAY RAL7038 1650 mm - Ø 10,5 mm		PC666








BIENAIR MICROMOTORS

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR MC3 HSH			 Spray regulation Swivel EXTENSION STRING LIGHT GRAY RAL7035 1500 mm - Ø 11 mm		PC666TM
MICROMOTOR MC3			 Spray regulation Swivel Safety Handpiece EXTENSION STRING ITAL GRAY RAL7038 1650 mm - Ø 10,5 mm		PC666SH
MICROMOTOR MX2			 Spray regulation Swivel EXTENSION STRING ITAL GRAY RAL7038 1500 mm - Ø 9,5 mm		PC734BXG




KAVO MICROMOTORS

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGTH / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR K192					PC124PKN
MICROMOTOR K190 K196 K201					PC135PKN
MICROMOTOR K700 K701 K703					PC148PKN




NSK MICROMOTORS

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR NLX NANO - PLUS ADEC AE53-NLZ			 Spray regulation Swivel LIGHT GRAY RAL7035 1800 mm - Ø 8,5 mm		PC601N
MICROMOTOR TIM 40 EL400			 Spray regulation Swivel EXTENSION STRING LIGHT GRAY RAL7035 1500 mm - Ø 9,5 mm		PC602N
MICROMOTOR M40 HSH			 Spray regulation Swivel EXTENSION STRING LIGHT GRAY RAL7035 1600 mm - Ø 11 mm		PC666TM



NSK MICROMOTORS

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR M40 HSH			 Spray regulation Swivel Safety Handpiece EXTENSION STRING ITAL GRAY RAL7038 1700 mm - Ø 11 mm		PC901NSH
MICROMOTOR NLX-NL400			 Spray regulation Swivel LIGHT GRAY RAL7035 1600 mm - Ø 10,5 mm		PC603N



FONA MICROMOTORS

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR BL					PC233F





PLANMECA MICROMOTORS

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR MINETTO					PC888P

FARO MICROMOTORS

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR MODO					PC756

TEKNEDENTAL MICROMOTORS

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR TKD DEFINITIVE LED			 Spray regulation Swivel EXTENSION STRING ITAL GRAY RAL7038 1650 mm - Ø 11 mm		PC668
MICROMOTOR TKD DEFINITIVE LED			 Spray regulation Swivel Safety Handpiece EXTENSION STRING ITAL GRAY RAL7038 1650 mm - Ø 11 mm		PC668SH