

artic digital

Carta de Formas



Artic® Digital

Salud oral en buenas manos



KULZER
MITSUI CHEMICALS GROUP

artic digital 6

Anteriores

Las siguientes formas estarán disponibles próximamente:

Anteriores superiores: O10, T15.

Anteriores inferiores: L05

O 20

46,1

9,7

8,4

S 35

50,6

10,8

9,2

O 30

46,6

10,9

8,4

S 40

52,3

12,4

9,8

R 20

42,9

9,9

7,8

T 10

44,1

9,0

8,2

T 20

45,6

9,8

8,8

R 40

51,1

13,3

9,5

T 30

45,0

10,8

8,1

S 05

40,3

8,9

7,4

T 35

47,4

10,0

8,8

S 10

44,2

9,1

7,9

T 40

49,1

12,1

8,7

S 20

44,0

9,9

7,9

S 30

48,2

10,9

8,6

artic digital 6

Anteriores

TO 20



43,6

10,2

7,7

TO 30



45,9

10,4

8,6

TO 40



47,8

11,0

8,9

L 10



33,5

7,8

5,1

L 15



33,9

9,5

5,2

L 20



35,2

9,1

5,4

L 30



35,9

9,0

5,5

L 40



36,0

9,0

5,5

L 45



38,9

9,6

5,9

L 50



39,2

11,1

5,9

L 55






40,2

10,6


6,0

Posteriores

De la angulación 20°, las siguientes tablillas estarás disponibles próximamente: 30U, 30L.

30 U	0°			27,2
30 L				7,9
32 U	0°			33,1
32 L				8,8
32 U	20°			31,1
32 L				9,1
34 U	20°			32,9
34 L				9,4
36 U	20°			34,7
36 L				10,0
				37,4
				9,0

	Anteriores superiores	Anteriores inferiores	Posteriores	0°	20°
	 O 20	L 30		30, 32	32
	 O 30	L 40		32	32
	 R 20	L 20		32	32, 34
	 R 30	L 25		32	36
	 R 40	L 55		32	36
	 S 05	L 10		30	32
	 S 10	L 10		30	32
	 S 20	L 20		30, 32	32
	 S 30	L 40		32	32, 34
	 S 35	L 55		32	36
	 S 40	L 55		32	36
	 T 10	L 10		30	32
	 T 20	L 30		30, 32	32
	 T 30	L 20		32	32
	 T 35	L 40		32	32
	 T 40	L 50		32	34
	 TO 05	L 10		30	32
	 TO 10	L 15		30	32
	 TO 20	L 20		30, 32	32
	 TO 30	L 40		32	32
	 TO 40	L 45		32	34



Esta carta de formas es conforme con la EN ISO 22112:2017.



Fabricante

Kulzer GmbH
Leipziger Straße 2
63450 Hanau
Alemania

Contacto en España

Kulzer S.A.
Marie Curie 19 bajo, oficina 5 ·
28521 Rivas Vaciamadrid, Madrid
infodental-es@kulzer-dental.com