



PIGRECO


PIGRECO è una lampada a luce diretta e indiretta con LED di ultima generazione a spettro solare. Grazie all'utilizzo di queste nuove sorgenti luminose, questa lampada riproduce fedelmente la qualità della luce naturale del Sole, rendendo ancora più elevato l'indice di resa cromatica.

Il suo design moderno e accattivante rende l'ambiente unico e originale, allo stesso tempo la sua luce di qualità che permette al dentista di operare in totale sicurezza e benessere.

PIGRECO is a direct and indirect light lamp with latest generation solar spectrum LED sources. Thanks to the use of these new light sources, this lamp faithfully reproduces the quality of natural sunlight, making the colour rendering index even higher.

Its modern and captivating design makes the environment unique and original, at the same time its quality light allows the dentist to operate in total safety and well-being.





Curva fotometrica
Photometric diagram

NEW

Nuovi LED a spettro solare che riproducono in modo identico la luce naturale, per una presa colore ancor più fedele!

New solar spectrum LEDs that reproduce natural light in an identical way, for an even more faithful colour capture!

Lampada fornita completa di kit di sospensione composto da cavi di sospensione metallici (pz. 4), cavo elettrificato trasparente (2000mm) e rosone tondo elettrificato (Ø100mm).

Lamp furnished with suspension kit made up of steel suspension cables (pcs. 4), transparent power cable (2000mm) and electrified round ceiling rose (Ø100mm).



Modelli / Models	PIGRECO	PIGRECO DIM
Codice / Code	364300074	364300075
Luce diretta* / Direct light*	100W 10508lm	
Luce indiretta* / Indirect light*	60W 6305lm	
Cablaggio / Wiring	ON/OFF	DIM DALI
Sorgente luminosa Luminous source	LED	
Temperatura di colore Colour temperature	5000K	
Indice di resa cromatica (IRC) Colour Rendering Index (CRI)	>97 - SUNLIGHT SPECTRUM	
Grado di protezione Protection degree	IP40	
Diffusore Diffuser	Metacrilato opale MOT Opal methacrylate MOT	
Peso / Weight	12 Kg	
Finitura / Finishing	Bianco opaco (RAL9003) / Matt white (RAL9003)	

