

PLANMECA DENTAL HOSES

The dental hoses for Tecnomed Italia dental units are made of silicone; they are Latex Free and comply with the UNI CEI EN ISO 13485 quality system.

The chemically stable gold-plated contacts ensure high electrical conductivity and maximum handpiece performance.

The supply tubes comply with REGULATION (EU) No. 10/2011 concerning plastic materials objects intended to meet food products and the exhaust tubes are resistant to the lubrication oils.

Tecnomed Italia dental hoses meet the strictest medical standards, complying with the essential requirements of the European Directives (93/42 / EEC) 2007/47 / EC on medical devices and 2000/95 / EC RoHs on the restriction of the use of certain specific substances in electrical and electronic equipment.

A complete range of dental hoses:

- Turbines
- Micromotors
- Scalers
- Syringes
- Freely Adaptable Dental Hoses

Features

- ⊕ High quality silicone sheath.
- ⊕ Latex free.
- ⊕ 12 months warranty

Technical data

- ⊕ Maximum pressure supply/return air 4,5 bar / 65 psi
- ⊕ Maximum pressure air/water spray 4 bar / 58 psi
- ⊕ Tensile strenght 70N



MICROMOTORS TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR MC3			RES. 83 Ω Spray regulation Swivel LIGHT GRAY RAL7035 890 mm - $\varnothing$ 10 mm 980 mm		PC812P
MICROMOTOR MC3			RES. 83 Ω Spray regulation Swivel LIGHT GRAY RAL7035 925 mm - $\varnothing$ 10 mm 1015 mm		PC812PL
MICROMOTOR MC3			RES. 825 Ω Spray regulation Swivel LIGHT GRAY RAL7035 890 mm - $\varnothing$ 10 mm 980 mm		PC812P2



MICROMOTORS TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGTH / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR MC3			RES. 825 Ω Spray regulation Swivel LIGHT GRAY RAL7035 925 mm - $\varnothing$ 10 mm 1015 mm		PC812PL2
MICROMOTOR K190 K196 K200			RES. 267 Ω Spray regulation Swivel LIGHT GRAY RAL7035 925 mm - $\varnothing$ 10 mm 1015 mm		PC823PL
MICROMOTOR K190 K196 K200			RES. 267 Ω Spray regulation Swivel LIGHT GRAY RAL7035 900 mm - $\varnothing$ 10 mm 1000 mm		PC823P

MICROMOTORS TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR K190 K196 K200			 Spray regulation Swivel LIGHT GRAY RAL7035 900 mm - Ø 10 mm 1000 mm		PC823P/2
MICROMOTOR MINETTO			 Spray regulation Swivel EXTENSION STRING LIGHT GRAY RAL7035 1800 mm - Ø 10 mm https://vimeo.com/480812578		PC888P

HIGHSPEED TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
HIGHSPEED HSH			RES. 825 Ω Spray regulation Swivel LIGHT GRAY RAL7035 925 mm - $\varnothing$ 10 mm 1015 mm		PC801P
HIGHSPEED HSH			RES. 430 Ω Spray regulation Swivel LIGHT GRAY RAL7035 825 mm - $\varnothing$ 10 mm 915 mm		PC802P
HIGHSPEED HSH			RES. 430 Ω Spray regulation Swivel LIGHT GRAY RAL7035 840 mm - $\varnothing$ 10 mm 930 mm		PC802PL


HIGHSPEED TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
HIGHSPEED HSH			RES. 430 Ω Spray regulation Swivel LIGHT GRAY RAL7035 925 mm - $\varnothing$ 10 mm 1015 mm	 	PC803PVL
HIGHSPEED HSH			RES. 560 Ω Spray regulation Swivel LIGHT GRAY RAL7035 915 mm - $\varnothing$ 10 mm 825 mm	 	PC803P

SYRINGE TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SYRINGE LUZZANI 3F			RES. 1,78 KΩ LIGHT GRAY RAL7035 910 mm - Ø 10 mm 1155 mm		PC845P

SCALER TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGTH / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SCALER EMS			RES. 1,33 KΩ Spray regulation Swivel LIGHT GRAY RAL7035 930 mm - $\varnothing$ 7 mm 1020 mm		PC834P
SCALER EMS O. F.			RES. 1,33 KΩ Spray regulation Swivel LIGHT GRAY RAL7035 930 mm - $\varnothing$ 8 mm 1020 mm		PC834PL
SCALER SATELEC O. F.			RES. 1,2 KΩ Spray regulation Swivel LIGHT GRAY RAL7035 930 mm - $\varnothing$ 8 mm 1020 mm		PC834PSL

SCALER TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGTH / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SCALER SATELEC			RES. 1,2 KΩ Spray regulation Swivel LIGHT GRAY RAL7035 930 mm - Ø 7 mm 1020 mm		PC834PS