

# OPEN END SCALER DENTAL HOSES

The dental hoses for Tecnomed Italia dental units are made of silicone; they are Latex Free and comply with the UNI CEI EN ISO 13485 quality system.

The chemically stable gold-plated contacts ensure high electrical conductivity and maximum handpiece performance.

The supply tubes comply with REGULATION (EU) No. 10/2011 concerning plastic materials objects intended to meet food products and the exhaust tubes are resistant to the lubrication oils.

Tecnomed Italia dental hoses meet the strictest medical standards, complying with the essential requirements of the European Directives (93/42 / EEC) 2007/47 / EC on medical devices and 2000/95 / EC RoHs on the restriction of the use of certain specific substances in electrical and electronic equipment.

A complete range of dental hoses:

- Turbines
- Micromotors
- Scalers
- Syringes
- Freely Adaptable Dental Hoses

- **High quality silicone sheath**
- **Latex free**
- **12 months warranty**

## **Technical data:**

Maximum pressure supply/return air: 4.5 bar / 65 psi

Maximum pressure air/water spray: 4 bar / 58 psi

Tensile strength: 70 N










# OPEN END CORDONS SCALER

HANDPIECE TYPE	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
EMS®		<p>EXTENSION STRING</p> <p>LIGHT GREY RAL7035</p> <p>L 2000 mm - Ø 7 mm</p>		PC314
EMS®		<p>EXTENSION STRING</p> <p>LIGHT GREY RAL7035</p> <p>L 2000 mm - Ø 7 mm</p>		PC314/B
EMS® O.F.		<p>EXTENSION STRING</p> <p>LIGHT GREY RAL7035</p> <p>L 2000 mm - Ø 8 mm</p>		PC314L




# OPEN END CORDONS SCALER

HANDPIECE TYPE	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
EMS® O.F.		<p>LIGHT GREY RAL7035</p> <p>L 1980 mm - Ø 8 mm</p>		PC314LW
EMS® I-PIEZON WHITE COLOR		<p>LIGHT GREY RAL7035</p> <p>L 1850 mm - Ø 8 mm</p>		PC314IPL
SATELEC® NSK VARIOS		<p>LIGHT GREY RAL7035</p> <p>L 2000 mm - Ø 7 mm</p>		PS115

# OPEN END CORDONS SCALER

HANDPIECE TYPE	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SATELEC® 3 WIRES		 <p>EXTENSION STRING</p> <p>WHITE RAL9010</p> <p>L 2000 mm - Ø 8 mm</p>		PS113
SATELEC® O.F. ULTRADENT NEWTRON		 <p>EXTENSION STRING</p> <p>LIGHT GREY RAL7035</p> <p>L 2000 mm - Ø 8 mm</p>		PS115L
NSK® VARIOS LUX		 <p>EXTENSION STRING</p> <p>LIGHT GREY RAL7035</p> <p>L 1850 mm - Ø 7 mm</p>	 <p>with LED 3.5 VDC</p>	PC118N

# OPEN END CORDONS SCALER

HANDPIECE TYPE	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MECTRON®		 <p>EXTENSION STRING</p> <p>LIGHT GREY RAL7035</p> <p>L 1700 mm - Ø 8 mm</p>		PS117MN