

DENTAL HOSES

The dental hoses for Tecnomed Italia dental units are made of silicone; they are Latex Free and comply with the UNI CEI EN ISO 13485 quality system.

The chemically stable gold-plated contacts ensure high electrical conductivity and maximum handpiece performance.







The supply tubes comply with REGULATION (EU) No. 10/2011 concerning plastic materials objects intended to meet food products and the exhaust tubes are resistant to the lubrication oils.







Tecnomed Italia dental hoses meet the strictest medical standards, complying with the essential requirements of the European Directives (93/42 / EEC) 2007/47 / EC on medical devices and 2000/95 / EC RoHs on the restriction of the use of certain specific substances in electrical and electronic equipment.







A complete range of dental hoses:



- Turbines
- Micromotors
- Scalers
- Syringes
- Freely Adaptable Dental Hoses

- ⊕ **High quality silicone sheath.**
- ⊕ **Latex free.**
- ⊕ **Medical device compliant with Directive 2007 / 47 / CE**
- ⊕ **Made in Europe Safety and quality**
- ⊕ **12 months warranty**

HANDPIECE TYPE	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
EMS®		 EXTENSION STRING LIGHT GREY RAL7035 L 2000 mm - Ø 7 mm		PC314
EMS®		 EXTENSION STRING LIGHT GREY RAL7035 L 2000 mm - Ø 7 mm		PC314/B
EMS® O.F.		 EXTENSION STRING LIGHT GREY RAL7035 L 2000 mm - Ø 8 mm		PC314L

HANDPIECE TYPE	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
EMS® I-PIEZON WHITE COLOR		 EXTENSION STRING LIGHT GREY RAL7035 L 1800 mm - Ø 8 mm		PC314IPL
EMS® I-PIEZON WHITE COLOR		 EXTENSION STRING LIGHT GREY RAL7035 Spray regulation L 1800 mm - Ø 8 mm		PC314IPLW
SATELEC® NSK VARIOS		 EXTENSION STRING LIGHT GREY RAL7035 L 2000 mm - Ø 7 mm		PS115

HANDPIECE TYPE	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SATELEC® 3 WIRES				PS113
SATELEC® O.F. ULTRADENT NEWTRON				PS115L
NSK® VARIOS LUX			 with LED	PC118N

HANDPIECE TYPE	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MECTRON®		 EXTENSION STRING LIGHT GREY RAL7035 L 1800 mm - Ø 8 mm		PS117MN