

OPEN END TURBINE DENTAL HOSES

The dental hoses for Tecnomed Italia dental units are made of silicone; they are Latex Free and comply with the UNI CEI EN ISO 13485 quality system.

The chemically stable gold-plated contacts ensure high electrical conductivity and maximum handpiece performance.

The supply tubes comply with REGULATION (EU) No. 10/2011 concerning plastic materials objects intended to meet food products and the exhaust tubes are resistant to the lubrication oils.

Tecnomed Italia dental hoses meet the strictest medical standards, complying with the essential requirements of the European Directives (93/42 / EEC) 2007/47 / EC on medical devices and 2000/95 / EC RoHs on the restriction of the use of certain specific substances in electrical and electronic equipment.

A complete range of dental hoses:

- Turbines
- Micromotors
- Scalers
- Syringes
- Freely Adaptable Dental Hoses

Features

- ⊕ High quality silicone sheath.
- ⊕ Latex free.
- ⊕ 12 months warranty

Technical data

- ⊕ Maximum pressure supply/return air 4,5 bar / 65 psi
- ⊕ Maximum pressure air/water spray 4 bar / 58 psi
- ⊕ Tensile strenght 70N

OPEN END CORDONS TURBINE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
HIGHSPEED LABORATORY HSH B2			<p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>250/1000 mm</p>		PL042B
HIGHSPEED HSH B2/3			<p>Spray regulation Swivel</p> <p>EXTENSION STRING ITAL GRAY RAL7038</p> <p>1650 mm - Ø 10,5 mm</p>		PC334
HIGHSPEED HSH B2/3			<p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>200/1000 mm</p>		PC534







OPEN END CORDONS TURBINE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
HIGHSPEED HSH M4			<p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1500 mm - Ø 11 mm</p>		PC444
HIGHSPEED HSH M4			<p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>200/1000 mm</p>		PC544
HIGHSPEED HSH M4 O.F.			<p>Spray regulation Swivel</p> <p>EXTENSION STRING</p> <p>LIGHT GRAY RAL7035</p> <p>1500 mm - Ø 11 mm</p>		PC655

OPEN END CORDONS TURBINE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
HIGHSPEED HSH M4 O.F.			<p>Spray regulation Swivel</p> <p>EXTENSION STRING</p> <p>ITAL GRAY RAL7038</p> <p>1650 mm - Ø 10,5 mm</p>		PC655TK
HIGHSPEED HSH CASTELLINI OLD			<p>Spray regulation Swivel</p> <p>EXTENSION STRING</p> <p>WHITE RAL9010</p> <p>1650 mm - Ø 10,5 mm</p>		PC306C
HIGHSPEED HSH CASTELLINI OLD			<p>Spray regulation Swivel</p> <p>EXTENSION STRING</p> <p>ITAL GRAY RAL7038</p> <p>1650 mm - Ø 10,5 mm</p>		PC306C/2

OPEN END CORDONS TURBINE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
HIGHSPEED HSH M4 O.F.			 <p>Spray regulation Swivel Safety Handpiece</p> <p>EXTENSION STRING</p> <p>ITAL GRAY RAL7038</p> <p>1650 mm - Ø 10,5 mm</p>		PC655TKSH
HIGHSPEED HSH M4 O.F. MC3 M40			 <p>Spray regulation Swivel</p> <p>EXTENSION STRING</p> <p>LIGHT GRAY RAL7035</p> <p>1500 mm - Ø 11 mm</p>		PC666TM