

FICHA TÉCNICA



CÁPSULAS DE AMALGAMA DENTAL PROCLINIC

REFERENCIAS	32120: 400 mg, 50 cápsulas 32121: 600 mg, 50 cápsulas
FABRICANTE / RESPONSABLE	CHL MEDICAL SOLUTIONS, Srl Via Felice Cavallotti, 15 20122 Milán - Italia
CLASIFICACIÓN	Producto Sanitario Clase IIa, según la Directiva 93/42/CEE relativa a productos sanitarios.
DESCRIPCIÓN	Cápsulas conteniendo aleación de metales y mercurio que, al mezclarlos, se usa la mezcla para sellar cavidades dentales.
COMPOSICIÓN	Cápsula de polietileno de baja densidad. Contenido: 2 partes: - 50% de Mercurio. - 50% de aleación de metales: plata, estaño, cobre, y puede contener trazas de zinc.
ESPECIFICACIONES	Cumple los requerimientos de ISO 24234 e ISO 13897. No Gamma 2 400 mg: cápsula de capuchón amarillo. 600 mg: cápsula de capuchón verde. Resistencia a la compresión (1 hora): >180MPa Resistencia a la compresión (24 horas): >460MPa Cambio dimensional: +0,08%
ENVASE	Envase conteniendo 50 cápsulas.
FINALIDAD PREVISTA	Para uso profesional dental. Las amalgamas dentales se utilizan para sellar las cavidades causadas por la caries dental.
VIDA ÚTIL	7 años.
ALMACENAMIENTO Y TRANSPORTE	Almacenar en lugar fresco y seco, por debajo de 25 °C. Proteger de la humedad.

TECHNICAL DATA SHEET



PROCLINIC DENTAL AMALGAM CAPSULES

REFERENCES	32120: 400 mg, 50 capsules 32121: 600 mg, 50 capsules
MANUFACTURER	CHL MEDICAL SOLUTIONS, Srl Via Felice Cavallotti, 15 20122 Milano - Italy
CLASSIFICATION	Medical Device class IIa, according to Directive 93/42/EEC concerning medical devices.
DESCRIPTION	Capsules containing alloy and mercury that after trituration, the mixture is used to seal the dental cavities.
COMPOSITON	Capsule shell: low density polyethylene. Contents: 2 parts: - 50% mercury - 50%: alloy of silver, tin, copper, and may contain small amounts of zinc.
SPECIFICATIONS	In conformity with the standard ISO 24234 and ISO 13897. Non Gamma 2. 400 mg: capsule with yellow cap. 600 mg: capsule with green cap. Compressive strength (1 hour): >180MPa Compressive strength (24 hours): >460MPa Dimensional change: +0,08%
PACKAGING	Package with 50 capsules.
INTENDED USE	For dental professional use only. Dental amalgams are intended to seal the cavities caused by tooth decay.
SHELF LIFE	7 years
HANDLING AND STORAGE CONDITIONS	Store in a cool dry place, below 25 °C. Protect from humidity.