

FICHA TÉCNICA

GUANTES LÁTEX SIN POLVO BLANCOS PROCLINIC		
REFERENCIAS	55010: Tamaño L 55020: Tamaño M 55030: Tamaño S 55040: Tamaño XS	
FABRICANTE / RESPONSABLE	CHL MEDICAL SOLUTIONS, Srl Via Felice Cavallotti, 15 20122 Milán - Italia	
CLASIFICACIÓN	Producto Sanitario Clase I, según la Directiva 93/42/CEE	
DESCRIPCIÓN	Guantes de látex sin polvo, no estériles, de un solo uso. Ambidiestros.	
COMPOSICIÓN	100% goma natural.	
ESPECIFICACIONES	Tamaño XS ≤ 80 mm Tamaño S 80 ± 10 mm Tamaño M 95 ± 10 mm Tamaño L 110 ± 10 mm	
ENVASE	Caja con 100 unidades	
VIDA ÚTIL	5 años	
FINALIDAD PREVISTA	Protección de barrera de la mano del usuario para la realización de exámenes médicos o clínicos, procedimientos diagnósticos y terapéuticos, en odontología y también para fines de laboratorio.	
INFORMACIÓN SOBRE LA SEGURIDAD	Los guantes están fabricados con látex, que puede causar reacciones alérgicas en algunas personas. Si esto ocurre, discontinuar su uso inmediatamente y consultar con el médico.	
ALMACENAMIENTO Y TRANSPORTE	Almacenar en lugar fresco y seco. Mantener alejado de la luz solar. Proteger de la humedad y evitar el calor excesivo.	
NORMATIVAS QUE CUMPLEN	EN 455:-1,-2,-3,-4 EN 420+A1 EN 374:-1	

TECHNICAL DATA SHEET

WHITE POWDER FREE LATEX GLOVES PROCLINIC		
REFERENCES	55010: Size L 55020: Size M 55030: Size S 55040: Size XS	
MANUFACTURER	CHL MEDICAL SOLUTIONS, Srl Via Felice Cavallotti, 15 20122 Milano - Italy	
CLASSIFICATION	Medical Device, Class I, according to Directive 93/42/EEC	
DESCRIPTION	Powder free latex gloves. Non-sterile. Ambidextrous. For single use only.	
COMPOSITION	100% natural rubber.	
SPECIFICATIONS	Size XS ≤ 80 mm Size S 80 ± 10 mm Size M 95 ± 10 mm Size L 110 ± 10 mm	
PACKAGING	100 units	
SHELF LIFE	5 years	
INTENDED USE	Gloves intended for use in the medical field to protect patient and user from cross contamination. Intended for conducting medical examinations, diagnostic and therapeutic procedures, in dentistry and in handling contaminated medical material.	
SAFETY INFORMATION	The gloves are made of natural rubber latex that may cause allergic reactions in some people. If allergic reactions occur, discontinue use immediately and consult the doctor.	
HANDLING AND STORAGE	Store in a cool dry place, away from direct sunlight. Protect from humidity and excessive heat.	
STANDARDS APPLIED	EN 455:-1,-2,-3,-4 EN 420+A1 EN 374:-1	