



FICHA TÉCNICA

GUANTES DE NITRIL PROCLINIC		
REFERENCIAS	7837: Tamaño L 7835: Tamaño M 7834: Tamaño S 7838: Tamaño XS	
FABRICANTE / RESPONSABLE	CHL MEDICAL SOLUTIONS, Srl Via Felice Cavallotti, 15 20122 Milán - Italia	
CLASIFICACIÓN	Producto Sanitario Clase I, según la Directiva 93/42/CEE	
DESCRIPCIÓN	Guantes de nitrilo, sin polvo, no estériles, de un solo uso. Ambidiestros. Color azul, microtexturizados.	
COMPOSICIÓN	100% nitrilo.	
ESPECIFICACIONES	Tamaño XS ≤ 80 mm Tamaño S 80 ± 10 mm Tamaño M 95 ± 10 mm Tamaño L 110 ± 10 mm	
ENVASE	Caja con 100 unidades	
VIDA ÚTIL	5 años	
FINALIDAD PREVISTA	Protección de barrera de la mano del usuario para la realización de exámenes médicos o clínicos, procedimientos diagnósticos y terapéuticos, en odontología y también para fines de laboratorio.	
INFORMACIÓN SOBRE LA SEGURIDAD	Producto de un solo uso. No usar si el envase o el guante están dañados.	
ALMACENAMIENTO Y TRANSPORTE	Almacenar en lugar fresco y seco. Mantener alejado de la luz solar. Proteger de la humedad y evitar el calor excesivo.	
NORMATIVAS QUE CUMPLEN	EN 455:-1,-2,-3,-4 EN 420+A1 EN 374:-1	

TECHNICAL DATA SHEET

NITRILE GLOVES PROCLINIC		
REFERENCES	7837: Size L 7835: Size M 7834: Size S 7838: Size XS	
MANUFACTURER	CHL MEDICAL SOLUTIONS, Srl Via Felice Cavallotti, 15 20122 Milano - Italy	
CLASSIFICATION	Medical Device, Class I, according to Directive 93/42/EEC	
DESCRIPTION	Nitrile gloves powdered free. Non-sterile. Ambidextrous. Blue colour. Microtextured. For single use only.	
COMPOSITION	100% nitrile.	
SPECIFICACIONES	Size XS ≤ 80 mm Size S 80 ± 10 mm Size M 95 ± 10 mm Size L 110 ± 10 mm	
PACKAGING	100 units	
SHELF LIFE	5 years	
INTENDED USE	Gloves intended for use in the medical field to protect patient and user from cross contamination. Intended for conducting medical examinations, diagnostic and therapeutic procedures, in dentistry and in handling contaminated medical material.	
SAFETY INFORMATION	For single use only. Do not use if package or glove are damaged.	
HANDLING AND STORAGE	Store in a cool dry place, away from direct sunlight. Protect from humidity and excessive heat.	
STANDARDS APPLIED	EN 455:-1,-2,-3,-4 EN 420+A1 EN 374:-1	