


Proclinic Expert Total Etch Adhesive



Proclinic Expert Total Etch Adhesive is a fifth-generation adhesive for acid-based, single-component total etching for enamel and dentine.

SHEAR RESISTANCE

Full dentine hybridisation produces a high resistance to adhesion.




TOTAL SEALING

Total sealing of dentinal tubules, which minimizes post-operative sensitivity.

BisGMA FREE

Contains no BisGMA, doing away with any unpleasant effects this may cause.

FLUORIDE RELEASE



Dr Carlos Gargallo
Clínica Carlos Gargallo
Registered number 3947

"Proclinic Expert adhesives have become an essential part of my daily clinical routine. They are easy to apply, non-drip and leave no uncomfortable residue after polymerisation. Feedback has been very positive and there has been no post-operative sensitivity. Without a doubt, they form a perfect complement for the new Proclinic Expert composites"

Total Etch Adhesive, 5 ml	36599
Total Etch Adhesive , 2x5 ml pack	36598


Proclinic Expert Self-Etch Adhesive



Proclinic Expert Self-Etch Adhesive

is a seventh-generation, fast-action, single-component system for selective etching and self-etching on dentine and enamel.

FLUORIDE RELEASE




* Compared to the three best-selling products on the Spanish market. Tests conducted by manufacturer.


SHEAR RESISTANCE

A high adhesive resistance ensures restorations stay in place.

Universal Self-Etch on enamel*



Other brands on enamel*



Universal Self-Etch on dentine*



Other brands on dentine*



- A 7% content of specially treated hydrophilic and hydrophobic fillers ensures a finished surface free from aqueous residues.
- Contains neither BisGMA nor HEMA, thus avoiding any skin irritations.

- Very thin films produce seamless edges between the restoration and the tooth.
- Colour changes as it sets (from brownish yellow to clear), indicating when to proceed with the next step.